

## Polishing Head PE.IU & PE.RU



Kemet P262 with  
PEIU polishing head

Automatic Preparation can be carried out using the PE.RU and PE.IU polishing heads in conjunction with any of the Kemet P230, P260 and P262 base units. The polishing heads are built strictly in accordance with EEC standards.

### PE.RU

This automatic polishing head applies a central pressure to 3, 4 or 6 samples Ø 25 to 40mm at a time. Ideal for large volume applications requiring extreme flatness

#### Specification

- Central pressure
- 6 samples Ø 20 to 30mm
- 4 samples Ø 40mm
- Control panel with sensitive keys
- Pressure 0.5 to 15 daN
- Rotation speed 20 to 100 rpm
- Time: 10 seconds to 99 minutes
- Voltage: 230V, single phase

### PE.IU

This automatic polishing head applies individual pressure to 1-4 samples at a time with Ø 25-40mm

#### Specification

- Individual pressure
- 1 to 4 samples Ø 20 to 40mm
- Pressure 0.5 to 7 daN
- Control panel with sensitive keys
- Rotation speed: 20 to 100 rpm
- Time: 10 seconds to 99 minutes
- Voltage: 230V, single phase